

SureFit™ Seating System

The Gillette SureFit™ Seating System is primarily used by children who have cerebral palsy or other neuromuscular conditions. Each child will receive a seating system that is configured specifically for him or her, using measurements provided by a therapist or seating specialist.

The SureFit™ Seating System features prefabricated plastic modules and foam cushions enclosed in slipcovers. Using a special interface, the seating system can be used with a wheelchair, stroller or SureFit™ floor base with casters. It's designed to easily adjust as a child grows.

For more detailed information about the SureFit™ Seating System, see [Using and Caring for a Custom Gillette Seat.](#)

To schedule an appointment, call 651-290-8707.



Make An Appointment

[651-290-8707](tel:651-290-8707)

Refer a Patient

[651-325-2200](tel:651-325-2200)

This information is for educational purposes only. It is not intended to replace the advice of your health care providers. If you have any questions, talk with your doctor or others on your health care team.

If you are a Gillette patient with urgent questions or concerns, please contact Telehealth Nursing at [651-229-3890](tel:651-229-3890).