

EDR053-SP

La órtesis de pie (FO) y la órtesis de pie tipo UCBL se usan para apoyar, distribuir la presión y ayudar a mantener la alineación correcta del pie, especialmente del talón y el arco.

Calzado

Es muy importante que use zapatos que le queden bien cuando usa una FO o una UCBL. Si sus zapatos no le ajustan bien, las órtesis no podrán funcionar correctamente y es posible que no sean cómodas.

Asegúrese de que su pie se ajuste cómodamente en el zapato con la órtesis dentro. En algunos casos, puede ser necesario un zapato un poco más grande. Los zapatos deportivos que se pueden amarrar o abrochar son una buena opción. Las FO y UCBL por lo general reemplazan las plantillas y se dejan en los zapatos en todo momento.

This information is for educational purposes only. It is not intended to replace the advice of your health care providers. If you have any questions, talk with your doctor or others on your health care team.

If you are a Gillette patient with urgent questions or concerns, please contact Telehealth Nursing at 651-229-3890.